



Product description

Order number	86.1509
Product description	Tube, suitable for Coombs test, material: glass, (LxØ): 75 x 11.5 mm, transparent, 294 piece(s)/StackPack

Product characteristics

Graduation	no
Application	Coombs test
Suitable for	Coombs test
Base shape	round base

Size

Diameter	11.5 mm
Length of product	75 mm

Material & colours

Product material	glass
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2352
Type of smallest subpackaging	StackPack
Piece(s) / inner box	2352
Piece (s) / outer case	2352
Piece(s) / pallet	56448
Depth of case	300 mm
Width of case	355 mm
Height of case	335 mm
Case volume	0.0357 cbm
Weight of product	0.0038 kg
Weight of case	11.22 kg
EAN of inner box	4038917086783
EAN case	4038917072366

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

